

Product datasheet for **GA213454**

Mouse Cd300lb activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd300lb CRISPRa kit - CRISPR gene activation of mouse CD300 molecule like family member B
Symbol:	Cd300lb
UniProt ID:	Q3U497
Synonyms:	CLM-7; Clm7; Gm253; LMIR5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)