

Product datasheet for **GA213460**

Mouse Trim47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim47 CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 47
Symbol:	Trim47
Components:	GA213460G1, Trim47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213460G2, Trim47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213460G3, Trim47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8C0E3
Synonyms:	2210023F24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)