

Product datasheet for GA213465

Mouse Tnrc6c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnrc6c CRISPRa kit - CRISPR gene activation of mouse trinucleotide repeat containing 6C
Symbol:	Tnrc6c
Components:	GA213465G1, Tnrc6c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213465G2, Tnrc6c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213465G3, Tnrc6c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UHC0</u>
Synonyms:	9930033H14Rik; mKIAA1582
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)