


Product datasheet for **GA213487**

Mouse Plekhd1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekhd1 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family D (with coiled-coil domains) member 1
Symbol:	Plekhd1
Components:	GA213487G1, Plekhd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213487G2, Plekhd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213487G3, Plekhd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	B2RPU2
Synonyms:	3830431G21Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)