

Product datasheet for **GA213489**

Mouse Sip1l1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sip1l1 CRISPRa kit - CRISPR gene activation of mouse signal-induced proliferation-associated 1 like 1
Symbol:	Sip1l1
Components:	GA213489G1, Sip1l1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213489G2, Sip1l1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213489G3, Sip1l1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C0T5</u>
Synonyms:	4931426N11Rik; AW213287; mKIAA0440; Spar
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)