

Product datasheet for GA213498

Mouse Cipc activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Cipc CRISPRa kit - CRISPR gene activation of mouse CLOCK interacting protein, circadian |
| Symbol: | Cipc |
| UniProt ID: | Q8R0W1 |
| Synonyms: | 2310044G17Rik; AI842544; AI853744; Kiaa1737; mKIAA1737 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)