

Product datasheet for **GA213517**

Mouse Wdr60 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr60 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 60
Symbol:	Wdr60
Components:	GA213517G1, Wdr60 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213517G2, Wdr60 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213517G3, Wdr60 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C761</u>
Synonyms:	D430033N04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)