

Product datasheet for GA213547

Mouse Wdr41 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr41 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 41
Symbol:	Wdr41
Components:	GA213547G1, Wdr41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213547G2, Wdr41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213547G3, Wdr41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UDP0
Synonyms:	B830029I03Rik; MSTP048
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

