


Product datasheet for **GA213582**

Mouse Pip4p1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pip4p1 CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4,5-bisphosphate 4-phosphatase 1
Symbol:	Pip4p1
Components:	GA213582G1, Pip4p1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213582G2, Pip4p1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213582G3, Pip4p1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TWL2</u>
Synonyms:	BC039161
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)