


Product datasheet for **GA213588**

Mouse Khnyn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Khnyn CRISPRa kit - CRISPR gene activation of mouse KH and NYN domain containing
Symbol:	Khnyn
Components:	GA213588G1, Khnyn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213588G2, Khnyn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213588G3, Khnyn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80U38</u>
Synonyms:	6330555F21; 9130227C08Rik; AI605259; Kiaa0323; mKIAA0323
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)