

Product datasheet for **GA213597**

Mouse Fam167a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam167a CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 167, member A
Symbol:	Fam167a
Components:	GA213597G1, Fam167a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213597G2, Fam167a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213597G3, Fam167a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P1G6</u>
Synonyms:	A030013D21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)