

Product datasheet for **GA213611**

Mouse Ranbp3l activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Ranbp3l CRISPRa kit - CRISPR gene activation of mouse RAN binding protein 3-like |
| Symbol: | Ranbp3l |
| Components: | GA213611G1, Ranbp3l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213611G2, Ranbp3l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213611G3, Ranbp3l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q6PDH4</u> |
| Synonyms: | C130037N17Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)