

Product datasheet for GA213625

Mouse C030006K11Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C030006K11Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA C030006K11 gene
Symbol:	C030006K11Rik
Components:	GA213625G1, C030006K11Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213625G2, C030006K11Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213625G3, C030006K11Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VE95</u>
Synonyms:	1110033113Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

