

Product datasheet for **GA213626**

Mouse Arhgap39 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arhgap39 CRISPRa kit - CRISPR gene activation of mouse Rho GTPase activating protein 39
Symbol:	Arhgap39
Components:	GA213626G1, Arhgap39 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213626G2, Arhgap39 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213626G3, Arhgap39 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59281</u>
Synonyms:	9530053N22; AI843066; D15Wsu169e; Kiaa1688; Porf-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)