

Product datasheet for GA213634

Mouse Mpped1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mpped1 CRISPRa kit - CRISPR gene activation of mouse metallophosphoesterase domain containing 1
Symbol:	Mpped1
Components:	GA213634G1, Mpped1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA213634G2, Mpped1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA213634G3, Mpped1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q91ZG2</u>
Synonyms:	6230403M19; AW456965; D15Bwg0669e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

