


Product datasheet for **GA213641**

Mouse Zbed4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zbed4 CRISPRa kit - CRISPR gene activation of mouse zinc finger, BED type containing 4
Symbol:	Zbed4
Components:	GA213641G1, Zbed4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213641G2, Zbed4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213641G3, Zbed4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80WQ9</u>
Synonyms:	1700009F06Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)