

## Product datasheet for **GA213645**

### Mouse Mroh2b activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Mroh2b CRISPRa kit - CRISPR gene activation of mouse maestro heat-like repeat family member 2B  |
| Symbol:       | Mroh2b  |
| Components:   | GA213645G1, Mroh2b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213645G2, Mroh2b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213645G3, Mroh2b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q7M6Y6</a>  |
| Synonyms:     | 4930455B06Rik; Heatr7b2; SPIF   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online »](#)