

Product datasheet for GA213671

Mouse Tmem44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem44 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 44
Symbol:	Tmem44
Synonyms:	1700007N03Rik; 9330161C17Rik; B230220N21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)