

## Product datasheet for **GA213681**

### Mouse Dzip3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Dzip3 CRISPRa kit - CRISPR gene activation of mouse DAZ interacting protein 3, zinc finger
Symbol:	Dzip3
UniProt ID:	<u><a href="#">Q7TPV2</a></u>
Synonyms:	2A-HUB; 2310047C04Rik; 6430549P11Rik; A230104G20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)