

Product datasheet for **GA213709**

Mouse Zfp523 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp523 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 523
Symbol:	Zfp523
Components:	GA213709G1, Zfp523 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213709G2, Zfp523 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213709G3, Zfp523 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BMU0
Synonyms:	2810027J07; BC025615; Znf76
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)