


## Product datasheet for **GA213763**

### Mouse Rprd1a activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Rprd1a CRISPRa kit - CRISPR gene activation of mouse regulation of nuclear pre-mRNA domain containing 1A                                                                                                                                                                                                      |
| Symbol:       | Rprd1a                                                                                                                                                                                                                                                                                                        |
| Components:   | GA213763G1, Rprd1a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213763G2, Rprd1a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213763G3, Rprd1a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q8VDS4</u>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | C77387; mKIAA4077                                                                                                                                                                                                                                                                                             |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)