

Product datasheet for **GA213766**

Mouse **Pik3c3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3c3 CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol 3-kinase catalytic subunit type 3
Symbol:	Pik3c3
UniProt ID:	<u>Q6PF93</u>
Synonyms:	5330434F23Rik; Vps34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)