

Product datasheet for **GA213815**

Mouse Trpm3 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Trpm3 CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily M, member 3 |
| Symbol: | Trpm3 |
| Synonyms: | 6330504P12Rik; 9330180E14; AU018608; B930001P07Rik; LTRPC3; MLSN2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)