

Product datasheet for GA213824

Mouse Cyp2c70 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp2c70 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily c, polypeptide 70
Symbol:	Cyp2c70
Components:	GA213824G1, Cyp2c70 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA213824G2, Cyp2c70 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA213824G3, Cyp2c70 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q91W64</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

