

Product datasheet for **GA213829**

Mouse Cyp2c23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp2c23 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily c, polypeptide 23
Symbol:	Cyp2c23
Components:	GA213829G1, Cyp2c23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213829G2, Cyp2c23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213829G3, Cyp2c23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A0A1D5RLI1
Synonyms:	AW107714; b2b386Clo; BC034834
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)