


Product datasheet for **GA213892**

Mouse Kctd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kctd3 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 3
Symbol:	Kctd3
Components:	GA213892G1, Kctd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213892G2, Kctd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213892G3, Kctd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BFX3</u>
Synonyms:	4930438A20Rik; 9330185B06; E330032J19Rik; NY-REN-45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)