

Product datasheet for **GA213910**

Mouse Mpp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mpp4 CRISPRa kit - CRISPR gene activation of mouse membrane protein, palmitoylated 4 (MAGUK p55 subfamily member 4)
Symbol:	Mpp4
Components:	GA213910G1, Mpp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213910G2, Mpp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213910G3, Mpp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6P7F1
Synonyms:	Dlg6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)