

Product datasheet for GA213911

Mouse Ino80d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ino80d CRISPRa kit - CRISPR gene activation of mouse INO80 complex subunit D
Symbol:	Ino80d
Components:	GA213911G1, Ino80d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213911G2, Ino80d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213911G3, Ino80d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q66JY2
Synonyms:	7330405I11; A430093A21Rik; BC010584; OTTMUSG00000002594
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)