

Product datasheet for **GA213926**

Mouse Ppip5k2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppip5k2 CRISPRa kit - CRISPR gene activation of mouse diphosphoinositol pentakisphosphate kinase 2
Symbol:	Ppip5k2
Components:	GA213926G1, Ppip5k2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213926G2, Ppip5k2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213926G3, Ppip5k2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6ZQB6
Synonyms:	AW555814; Cfap160; D330021B20; Hisppd1; mKIAA0433; Vip2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)