

Product datasheet for **GA213939**

Mouse Tmem203 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem203 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 203
Symbol:	Tmem203
Components:	GA213939G1, Tmem203 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213939G2, Tmem203 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213939G3, Tmem203 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8R235
Synonyms:	C730025P13Rik; Hbebp1; RP23-132N23.14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)