


## Product datasheet for **GA213940**

### Mouse *Lrrc26* activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                                    |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                                       |
| Description:  | <i>Lrrc26</i> CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 26                                                                                                                                                                                                                                      |
| Symbol:       | <i>Lrrc26</i>                                                                                                                                                                                                                                                                                                                      |
| Components:   | GA213940G1, <i>Lrrc26</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213940G2, <i>Lrrc26</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213940G3, <i>Lrrc26</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q91W20</a>                                                                                                                                                                                                                                                                                                             |
| Synonyms:     | RP23-132N23.19                                                                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                            |

#### Product images:



[View online »](#)