

Product datasheet for **GA213955**

Mouse Ddx31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx31 CRISPRa kit - CRISPR gene activation of mouse DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 31
Symbol:	Ddx31
Components:	GA213955G1, Ddx31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213955G2, Ddx31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213955G3, Ddx31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NZQ2</u>
Synonyms:	5830444G11Rik; Gm997
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)