

Product datasheet for **GA213968**

Mouse Plpp7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plpp7 CRISPRa kit - CRISPR gene activation of mouse phospholipid phosphatase 7 (inactive)
Symbol:	Plpp7
UniProt ID:	Q91WB2
Synonyms:	D830019K17Rik; NET39; Ppapdc3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)