

Product datasheet for **GA213971**

Mouse Pip5k11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pip5k11 CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4-phosphate 5-kinase-like 1
Symbol:	Pip5k11
UniProt ID:	Q6U7H8
Synonyms:	BC028795
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)