

Product datasheet for **GA213983**

Mouse Pkp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pkp4 CRISPRa kit - CRISPR gene activation of mouse plakophilin 4
Symbol:	Pkp4
UniProt ID:	<u>Q68FH0</u>
Synonyms:	5031422I09Rik; 9430019K17Rik; Armrp; p0071
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)