

Product datasheet for **GA213989**

Mouse Atp5g3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Atp5g3 CRISPRa kit - CRISPR gene activation of mouse ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit C3 (subunit 9)
Symbol:	Atp5g3
Components:	GA213989G1, Atp5g3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213989G2, Atp5g3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213989G3, Atp5g3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56384</u>
Synonyms:	6030447M23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)