

Product datasheet for **GA213994**

Mouse P2rx3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2rx3 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel, 3
Symbol:	P2rx3
UniProt ID:	Q3UR32
Synonyms:	4930513E20Rik; P2X3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)