

Product datasheet for **GA214003**

Mouse Slc35c1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Slc35c1 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 35, member C1 |
| Symbol: | Slc35c1 |
| Components: | GA214003G1, Slc35c1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214003G2, Slc35c1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214003G3, Slc35c1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | Q8BLX4 |
| Synonyms: | E430007K15Rik; fuct1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)