

Product datasheet for GA214005

Mouse Prrg4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prrg4 CRISPRa kit - CRISPR gene activation of mouse proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane)
Symbol:	Prrg4
UniProt ID:	<u>Q8BGN6</u>
Synonyms:	9930111118Rik; TMG4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)