

## Product datasheet for GA214008

### Mouse Zfp770 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Zfp770 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 770                                                                                                                                                                                                                                  |
| <b>Symbol:</b>       | Zfp770                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA214008G1, Zfp770 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214008G2, Zfp770 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214008G3, Zfp770 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <a href="#">Q8BIQ8</a>                                                                                                                                                                                                                                                                                        |
| <b>Synonyms:</b>     | 6430601A2IRik; Znf770                                                                                                                                                                                                                                                                                         |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

