

Product datasheet for GA214019

Mouse Sptlc3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Sptlc3 CRISPRa kit - CRISPR gene activation of mouse serine palmitoyltransferase, long chain base subunit 3 |
| Symbol: | Sptlc3 |
| UniProt ID: | Q8BG54 |
| Synonyms: | C130053K05Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)