


Product datasheet for **GA214065**

Mouse Tm4sf4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tm4sf4 CRISPRa kit - CRISPR gene activation of mouse transmembrane 4 superfamily member 4
Symbol:	Tm4sf4
Components:	GA214065G1, Tm4sf4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214065G2, Tm4sf4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214065G3, Tm4sf4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	lItmp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)