


Product datasheet for **GA214067**

Mouse Clnr1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Clnr1 CRISPRa kit - CRISPR gene activation of mouse clarin 1
Symbol:	Clnr1
Components:	GA214067G1, Clnr1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214067G2, Clnr1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214067G3, Clnr1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8K445
Synonyms:	A130002D11Rik; clarin-1; USH3; Ush3a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)