


Product datasheet for **GA214115**

Mouse Plppr4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plppr4 CRISPRa kit - CRISPR gene activation of mouse phospholipid phosphatase related 4
Symbol:	Plppr4
Components:	GA214115G1, Plppr4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214115G2, Plppr4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214115G3, Plppr4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TME0</u>
Synonyms:	A330086D10; D3Bwg0562e; Lppr4; mKIAA0455; PRG-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)