

## Product datasheet for GA214121

### Mouse *Kyat3* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Kyat3 CRISPRa kit - CRISPR gene activation of mouse kynurenine aminotransferase 3
Symbol:	Kyat3
UniProt ID:	<a href="#">Q71RI9</a>
Synonyms:	Ccbl2; Kat3; KATIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)