


Product datasheet for **GA214130**

Mouse Phf24 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Phf24 CRISPRa kit - CRISPR gene activation of mouse PHD finger protein 24
Symbol:	Phf24
Components:	GA214130G1, Phf24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214130G2, Phf24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214130G3, Phf24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80TL4</u>
Synonyms:	GINIP; N28178
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)