

Product datasheet for GA214137

Mouse Slc25a51 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a51 CRISPRa kit – CRISPR gene activation of mouse solute carrier family 25, member 51
Symbol:	Slc25a51
Components:	GA214137G1, Slc25a51 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA214137G2, Slc25a51 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA214137G3, Slc25a51 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q5HZI9</u>
Synonyms:	9130208E07Rik; D130005A03Rik; Gm138; Mcartl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

