

## Product datasheet for **GA214142**

### Mouse Zfp189 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Zfp189 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 189                                                                                                                                                                                                                                  |
| Symbol:       | Zfp189                                                                                                                                                                                                                                                                                                        |
| Components:   | GA214142G1, Zfp189 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214142G2, Zfp189 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214142G3, Zfp189 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | BB131811; C430015I23Rik                                                                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)