

Product datasheet for GA214147

Mouse Ptbp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptbp3 CRISPRa kit - CRISPR gene activation of mouse polypyrimidine tract binding protein 3
Symbol:	Ptbp3
UniProt ID:	Q8BHD7
Synonyms:	5830471K22Rik; AA407443; AI462022; AW107884; C86549; Rod1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)