

Product datasheet for **GA214173**

Mouse *Lrrc41* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Lrrc41</i> CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 41
Symbol:	<i>Lrrc41</i>
Components:	GA214173G1, <i>Lrrc41</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214173G2, <i>Lrrc41</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214173G3, <i>Lrrc41</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8K1C9
Synonyms:	AA409966; AW555107; D630045E04Rik; D730026A16Rik; MUF1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)